



BOSTON XO

Boston XO is a tough, stable, high Dk material that can be made into a wide variety of special designs for demanding visual needs. Boston XO has an ISO/Fatt Dk of 100 with a proven record of international success as a premium extra oxygen exceptional stability GP material.

Indications/Features

- Myopia
- Hyperopia
- Astigmatism
- Presbyopia
- Daily wear



Material	hexafon A	Saline Absorption	1.0 %
Oxygen Permeability Dk	100*	Refractive Index	1.415
Wetting Angle (Captive Bubble)	49 deg.	Specific Gravity	1.27
Hardness (Rockwell R)	112	Tint	Blue

