



STARSOFT 75 UV

Contaflex's high-water content means greater permeability, leading to improved user comfort and excellent clarity.

The water content in Contaflex 75% is ideally suited to bandage and other therapeutic lens applications.

Technical information

Material Characteristics:

Material	Contaflex 75 %	Light transmission	95%
Water Content	74.5%	CE Approved	Yes
UV Blocker	Yes	Use time	6 months
Handling tint	Clear	Material belonging according to the FDA	III
DK	43x10 ¹¹	Manufacturing	Lathing

